# TECHNISCHE INFORMATIONEN

# TECHNICAL INFORMATION

01-2025

Rodenstock Eyewear

Porsche Design

# INHALT

## CONTENTS

1.	<b>Fassungsmaterialien</b>	5
	Frame materials	
2.	<b>Pflegehinweise</b>	6
	Care instructions	
3.	<b>Kennzeichnung</b>	8
	Identifictaion	
4.	<b>Verglasungs- und Anpassungshinweise</b>	9
	Glazing and adjustment notes	
4.1	<b>Rodenstock T-Lite &amp; Porsche Design Ultra Light</b>	10
	Rodenstock T-Lite & Porsche Design Ultra Light	
4.2	<b>RXP® (Polyamid)</b>	10
	RXP® (Polyamide)	
4.3	<b>Kunststoffaugenrand</b>	11
	Plastic eye rim	
4.4	<b>Geteilte Backe</b>	14
	Splitted lug	
4.5	<b>Randlosfassungen mit Dübelsteckverbindung</b>	17
	Rimless frames with dowel plug-in connection	
4.6	<b>Rodenstock – One Strip of Metal</b>	20
	Rodenstock – One Strip of Metal	
4.7	<b>Büffelhorn</b>	22
	Buffalo horn	
5.	<b>Scharniere</b>	23
	Hinges	
5.1	<b>Rodenstock Roll-Hinge</b>	24
	Rodenstock Roll-Hinge	
5.2	<b>Titan-Zylinder-Scharnier</b>	25
	Titanium-Cylinder-Hinge	
5.3	<b>Rodenstock „Orbit“</b>	26
	Rodenstock "Orbit"	
5.4	<b>Rodenstock „Loop“</b>	27
	Rodenstock "Loop"	

# INHALT

## CONTENTS

5.5	<b>Porsche Design „Ball Tec“</b>	30
	Porsche Design "Ball Tec"	
5.6	<b>Porsche Design „Motion“</b>	32
	Porsche Design "Motion"	
5.7	<b>Porsche Design „Absorber“</b>	35
	Porsche Design "Absorber"	
6.	<b>Sonnenschutzgläser</b>	38
	Sun Protection Lenses	
6.1	<b>Funktion</b>	39
	Function	
6.2	<b>Gläser</b>	39
	Lenses	
6.3	<b>Materialien</b>	40
	Materials	
6.4	<b>Filterkategorien</b>	41
	Filter categories	
7.	<b>Tipps und Hinweise</b>	42
	Tips and hints	
7.1	<b>Bügelenden (Schrumpfschlauch)</b>	43
	Temple tips (Heat-shrink tubing)	
7.2	<b>Glasbefestigung mit Nieten</b>	44
	Lens fixation with rivets	
8.	<b>Service</b>	46
	Services	

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## Stainless Steel

Stainless steel is a very strong material and is characterized as a frame material mainly by its durability and robustness, as well as its very good flexibility and spring properties. A stainless steel spectacle frame can generally withstand high levels of stress.

## RXP® - Refined Excellence in Polyamide

The high-performance Polyamide RXP® is produced in Switzerland by high-tech injection moulding. This innovative frame material is very suitable to create very fine and delicate frames in a wide colour range. It is ultra lightweight, highly flexible but still consistent in form.

Additional features are: high impact and scratch resistant, easy to adjust and very skin friendly.

## Cellulose Acetate

One outstanding feature of this plastic is its excellent dermatological properties, which make it feel very comfy against the skin. Rodenstock mainly uses Italian acetate available in a range of impressive colour variants and finishes. Simply by gently heating and adjusting, acetate can very easily and effectively be tailored to create the right fit.

## Horn

Rodenstock uses selected horn from the Indian water buffalo, which is turned into an exceptionally high-quality, hypoallergenic frame material in an elaborate manufacturing process. Color variations in natural materials can occur and are a characteristic of the natural product. They make each piece unique.

## Gold

Gold is a highly prized precious metal that symbolizes enduring value. Thanks to their beautiful luster and rarity, coupled with their robustness, gold-framed spectacles are just the thing for those looking for something truly special. Every Rodenstock gold frame is a real masterpiece that involves more than 500 production steps carried out by hand.

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## Titanium

Rodenstock titanium frames are made of pure titanium or beta-titanium. These materials stand out thanks to their extremely impressive corrosion resistance and their high level of strength, which enables amazingly lightweight metal frames. Where the frame design creates special requirements in terms of flexibility, for example with small cross sections, beta-titanium is used instead of pure titanium: a very strong and completely nickel-free titanium alloy.

As what is known as a 'biocompatible' material, titanium has particularly good dermatological properties. Titanium frames are especially recommended for spectacle-wearers who are allergic to nickel. In this respect, Rodenstock titanium frames offer the highest levels of safety.

### 1. General

In order to enjoy your product as long as possible, please note the following care instructions:

- All Rodenstock spectacles can be cleaned with a standard microfiber cloth. Dry cleaning should be avoided if possible, because dust particles on the frame can be rubbed into it and cause scratches
- Spectacles should be rinsed under lukewarm running water with a pH-neutral cleaning agent, a diluted dishwashing detergent that has no relubricating properties, or a solvent-free spectacle care product. A clean, fine microfiber or cotton cloth is recommended for drying
- Do not use harsh cleaners, (alcohol-based) glass cleaners, liquids containing solvents (acetone, etc.), acids, alkalis or other aggressive chemical cleaners to clean the spectacles
- Protect the spectacles from high levels of moisture and temperatures above 55 °C, e.g. in the sauna, steam bath or in a car parked in the sun. These could have a negative effect on the spectacles

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## **2. Buffalo horn**

To ensure that the natural product lasts as long as possible, the following important care instructions should be considered:

- Avoid high temperatures above 55°C for long periods of time.
- Occasional applications of the provided care cream and subsequent polishing of the frame with a soft cloth protects against embrittlement.
- Excessive moisture, aggressive detergents or cleaning in an ultrasonic bath can cause irreparable damage.
- Please have your natural horn frame checked regularly by your optician.

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\*2 = Cat. 2

\*3 = Cat. 3

Every model of Rodenstock spectacles comes with a unique identifier to make it easier for you to work with our spectacles.

#### 1. Material identifier

a. Acetate or polyamide:

If the frame is made from either polyamide or acetate, you can check by looking at the left nose pad, where you will find a "CA" for acetate or a "RXP®" for polyamide.

b. Titanium:

If a frame is made from titanium, this is indicated on the right temple. On prescription frames this can also be discerned by checking the imprint of the demo lenses.

#### 2. Size indicator

By checking the ▲ or ▼ prior to the lens width you can tell at a glance whether there are smaller or larger sizes of this model.

▲ indicates that a larger size exists

▼ indicates that a smaller size exists

- indicates that there is only one size

On plastic frames the size identifier is indicated in the right nose pad; on metal/titanium frames it is positioned on the back of the bridge or on the left temple, as well as the temple length.

#### 3. Lens Filter identifier

All our sunglasses are marked with additional letters to identify the lens filter. You will find these letters, printed inside of the temple, directly behind the filter category.

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## 4.1 RODENSTOCK T-LITE & PORSCHE DESIGN ULTRA LIGHT

Glazing this highly flexible titanium frame is comparable to glazing a nylon thread frame. To achieve full stability and protect the lenses, we recommend:

- using plastic lenses with a refractive index of 1.60 or greater

- a minimum lens thickness of 2.2 mm

- Dimensions:

· T-Lite

thickness T-rim: 0.60 mm

groove depth: 0.50 mm

· Ultra Light

thickness rim wire: 0.90 mm

groove depth: 0.40 mm

- for T-lite frames:

· if necessary, widen or deepen the groove in the area of the T-profile slightly to ensure that the T-profile lies fully on the flat facet

· only lightly chamfer the outside edges in the area of the T-profile so that the largest possible bearing surface remains

- chamfering the groove edge after edging and grooving

- cleaning the lens with ultrasound after processing in order to completely remove all processing residue from the groove

For this model, it is not necessary to modify the frame due to increased lens radii because of the extremely flexible design of the rim.

## 4.2 RXP® (POLYAMIDE)

- For thin-walled polyamide fronts, we recommend using the original optician former to determine the lens shape, instead of tracing the rim

- Insert the lenses without heating the frame

- Avoid high temperatures or uneven heating as they may cause irreversible deformation

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Heat polyamide temples (without a metal insert) in the adjustable area uniformly up to approx. 140°C. Depending on the thickness of the temples, this takes 1–2 minutes. The correct temperature has been reached if the temples start to deform when they are tapped lightly. Then shape them into the desired shape and allow to cool. Adjustment temperatures above 160 °C can damage the material and must be avoided. Do not heat the temples near the joints and decorative elements, as this can lead to damage.

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### 4.3 PLASTIC EYE RIM

For this glazing a tweezer and a screwdriver are required.

Step 1:

- In the first step remove the rim lock screw. Turn the hinge screw counterclockwise and place it in a safe place
- Please be careful while removing the screws. Make sure that the screwdriver is not larger than the slot of the screw head

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
Step 2:

- For removing the eye rim including the demo lenses, the locking block screw must be opened. You have to turn the screw counterclockwise
- With the help of tweezers, the screws can be removed more easily
- Note that the two screws should not be interchanged.
- Remove the temple. It is important that the two washers remain on the temple

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#### Step 3:

- In the next step the plastic ring with the demo lens is removed from the spectacle frame. Then the demo lens can be pressed out of the plastic ring
- For removing the demo lens from the plastic eye rim, start at the center of the top of the lens and remove the ring. Be careful while pushing it from the inside to the outside

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#### Step 4:

- After the lens has been cut, the frame can now be reassembled. Make sure that both screws and the two washers are still present

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
#### Step 5:

- The lenses are inserted into the ring from behind without heating it. The glazing is started at the nose area
- Then proceed piece by piece with even pressure, so that the lens is centered in the plastic ring

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

#### Step 6:

- A marking guide is located in the rim lock area for precise positioning of the ring
- The screw can be threaded more easily with the help of the tweezers
- In order to prevent the screwdriver from slipping, it is useful to work on a firm underground. Turn the rim lock screw clockwise

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
#### Step 7:

- Do not close the rim lock completely to insert the temple
- Insert the temple incl. washers and fix it with the hinge screw
- Close the rim lock completely and adjust the hinge friction

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
## 4.4 SPLITTED LUG

For this glazing a tweezer and a screwdriver are required.

#### Step 1:

- In the first step remove the rim lock screw. Turn the screw counterclockwise and place it in a safe place.
- Please be careful while removing the screws. Make sure that the screwdriver is not larger than the slot of the screw head

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### Step 2:

- For removing the demo lens, the hinge screw must be opened. You have to turn the screw counterclockwise.
- With the help of tweezers the screws can be removed easily.
- Note that the two screws should not be interchanged.

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
### Step 3:

- Remove the temple from the guide. It is important that the two washers remain on the temple
- Then the demo lens can be pressed out of the spectacle frame

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
### Step 4:

- After the lens has been cut, the frame can now be reassembled. Make sure that both screws and the two washers are still present. With some oil the washers can be fixed.

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
Step 5:

- Insert the rim lock screw, but do not close it completely, to insert the temple
- Insert the temple incl. washers and fix it with the hinge screw
- Close the rim lock completely and adjust the hinge friction

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

## 4.5 RIMLESS FRAMES WITH DOWEL PLUG-IN CONNECTION

The technical notes refer to all rimless frames with dowel plug-in connection. Double-pin plug sleeves, dowel cutting pliers and assembly pliers are required for glazing.

The design defines the placement of the drill holes. The drill holes can be made either manually or with the automatic lens drill.

### Step 1:

According to the type of application, please refer to the illustration for the technical data required for this. The most important measurement units are visualized on this chart.

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


### Step 2:

- There are standardized distances for the positioning of the holes
- The following principle applies to all our rimless frames: The distance between the two drill holes is always 3.0 mm
- Distance between the outer hole and the edge of the lens is always 2.6 mm

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
### Step 3:

- Preparation of the lenses. When ordering the lenses ensure with your lens supplier that the lens thickness at the outside drill holes is minimum 2.2 mm and maximum 6.4 mm
- Cut the lenses to the exact shape and drill the holes with a diameter of 1.40 mm
- Drill each pair of parallel holes from the convex front side of the lenses at an angle of 90 degrees to the front surface of the lenses. Slightly chamfer the holes on the front and rear side

### Step 4:

- The frame parts are mounted safely and easily using standard high-quality tools for rimless frames. Insert the sleeves into the drill holes on the back of the lens

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
#### Step 5:

- With the help of a demo lens press against the traverse of the sleeves
- Cut off the sleeves perfectly flush with the front surface of the lens
- Pre-position the frame part at the front side of the lens

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
#### Step 6:

- Press the demo lens against the traverse of the sleeves to avoid any gap while pressing in by hand the first holding ring of the pins
- Position the pliers precisely onto the traverse of the sleeves and press against the traverse until it is flush with the lens. Keep this position under permanent pressure while closing the pliers and pressing the pins into the sleeves in one step

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
## 4.6 RODENSTOCK – ONE STRIP OF METAL

For this glazing a pull-in tape, scissors and flat pliers are required.

### Step 1:

Push the pull-in tape between the metal rim of the frame and the lens in the area where the nylon thread is attached. To make this easier, the end of the tape can be cut with scissors so that it forms a long, narrow point.

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



### Step 2:

- Hold the frame with the lens well
- Grab the tape with flat pliers and pull it halfway through

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### Step 3:



- Hold both ends of the pull-in tape and pull it down along the groove. Pull the nylon thread over the lens edge until the lens jumps out
- Please note that you only hold the frame so that the metal insert does not come off

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### Step 4:

- When glazing, make sure that the metal insert is positioned in the lens groove. Then thread the nylon thread into the groove using the pull-in tape

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## 4.7 BUFFALO HORN

### Glazing:

- Match the lens shape and size exactly to the frame shape.
- Lens groove depth for full rim frames is approx. 0.5 mm.
- Break the tip of the glass bevel slightly before insertion to avoid tension.
- Natural horn cannot be cold glazed.
- The lenses are inserted from the front.

### Heating:

- Before heating, lightly applying care cream or a little vegetable oil to the inner edge of the frame and the glass bevel can make insertion easier. This also protects the surface structure of the frame.
- Natural horn absorbs heat slowly and transfers it inwards very slowly. It is therefore advisable to apply heat gradually.
- Set the temperature on the valve to approx. 200 - 220 °C.
- Heat for 5-8 seconds, then interrupt for approx. 5 seconds.
- Allow the heat to soak in to soften the material.
- Repeat the process two to a maximum of three times until the horn becomes soft.

### Please note:

- Natural horn cannot be cold-fitted.
- Proceed as described in the "Heating" chapter.
- The titanium pad arms can be adapted to any anatomy.

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
## 5.1 RODENSTOCK ROLL HINGE

### Features of the Roll-Hinge:

- Smart solution for hinge axis
- Fully integrated RXP® hinge
- Extremely even and soft hinge motion
- Individually adjustable
- Durable and maintenance-free

### Temple Dismounting

- The Roll Hinge can be opened with standard tools
- Fix upper screw
- Open lower screw
- Push out upper screw with the lower screw

Schraubensicherung/  
Screw lock

### Temple Mounting

- Place the upper screw with the anti-turn device in the notch of the temple
- Fix in place the upper screw and center within the hinge
- Screw in the lower screw and adjust the hinge friction

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
## 5.1 TITANIUM-CYLINDER-HINGE

### Features of the Titanium-Cylinder-Hinge:

- Screwless hinge
- Extremely even hinge motion
- Self-adjusting
- Durable and maintenance-free

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**Do not use solvents or intense ultrasonic bath for cleaning the frames. Best for cleaning is water with mild cleaning agent.**

In the rare case of poor hinge function do not disassemble the hinge. Please follow the following steps:

- Hinge is too loose: Order a new metal front or a new lug-temple-combination
- Hinge is too tight: Try to reestablish the initial hinge friction by inserting a small drop of oil into the hinge. Then open and close it several times in order to distribute the oil evenly in the hinge

If the hinge function is still poor order a new metal front or a new lug-temple-combination.

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**Do not use solvents or intense ultrasonic bath for cleaning the frames. Best for cleaning is water with mild cleaning agent.**

In the rare case of poor hinge function do not disassemble the hinge. Please follow the following steps:

- Hinge is too loose: Order a new frame or a new lug-temple-combination
- Hinge is too tight: Try to reestablish the initial hinge friction by inserting a small drop of oil into the hinge. Then open and close it several times in order to distribute the oil evenly in the hinge

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## 5.4 RODENSTOCK LOOP

### Features of the Loop-Hinge:

- Smart hinge solution
- Even hinge motion
- Individually adjustable with standard tools

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### Temple mounting Loop-Hinge

To replace the temples, you will need tweezers, a Phillips screwdriver and a flat-head screwdriver. The flat-head screwdriver should be relatively wide and cover the entire screw head.

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#### Step 1:

In the first step, remove the threaded sleeve. Turn it anticlockwise and place the screw in a safe place. For easy handling, place the glasses on a firm surface.

Some of the screw can only be loosened with a lot of pressure. The screw can be removed more easily with the help of tweezers.

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Take care when removing the threaded sleeve. Make sure that the screwdriver is the same size as the slotted threaded sleeve.

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#### Step 2:

- To remove the screw sleeve, it must be pushed out from behind using a Phillips screwdriver or similar.
- Place it in a safe place and now remove the temple from the guide.
- Make sure that the washer disks do not get lost.

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
Step 3:

When mounting the temple, make sure that the square screw head is completely countersunk in its guide.

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If the screw cannot be pushed in completely, you can turn the screw head slightly until it disappears completely into the hole.

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Step 4:

Use the tweezers to insert the threaded sleeve with the slotted side facing downwards. Place the glasses on a stable surface so that you can tighten the threaded sleeve with a little pressure and prevent it from slipping. Turn the screw clockwise.

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
## 5.5 PORSCHE DESIGN BALL TEC

### Features of the Ball Tec hinge:

- Screw-less hinge
- Extremely even hinge motion
- Durable and maintenance-free
- Temples exchangeable by opticians


Please note:

Do not lubricate the hinge! Hinge function can be impaired by oil or other unsuitable lubricants.

### Dismounting of the temple

Step 1:

- Hold the front quite close to the hinge
- Hold the temple at the hinge like shown in the picture
- Press the temple against the lug until the ball snaps out

### Dismounting of the ball

Step 1:

- Open the temple slightly (approximately 45°). If the temple is already broken, the position does not matter, as only the ball is inside the front

Step 2:

- The tool must be placed on a flat surface. Place the frame on the tool and make sure that the surface of the ball is positioned on the pin of the tool. Hold the frame with both hands and press on the lug (next to the ball) with your thumbs on both sides

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**Step 3:**

- Then the ball or temple pops out. Please be careful not to damage the front by positioning the tool incorrectly. Furthermore, be careful and do not injure yourself by the ball popping out

Order number (Dismounting Tool Temple Ball Tec Series): MPOBALLT001


**Mounting of the temple**

**Step 1:**

- Align the flat parts of the ball in the direction of the slot in the lug
- Insert the ball into the slot as deep as it is possible without applying force

**Please note:**

No additional lubrication is required when exchanging the temple. Hinge function can be impaired by using oil or other unsuitable lubricants.

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Step 2:

- Place the lug with the prepositioned ball on a several times folded cleaning cloth (temple not in totally opened position)
- Firmly press against the ball and the temple while moving the temple slightly back and forth until the ball snaps into place

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### Dismounting of temple


Step 1:

- Move temple to closed position.
- Turn the hinge screw counterclockwise, remove it and store it in a safe place.  
(With the help of tweezers, the screw can be removed.)

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### Step 2:

- Move temple in closed position backwards from the lug until the temple is completely free.
- Check, if the hinge washer remains on the temple. If so, the hinge washer can be removed with a small screwdriver or tweezers and store it in a safe place.

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## Mounting of the temple

### Step 1:

- Place hinge washer on screw tube.
- Place one small drop glue (e.g. Loctite 242 Threadlocker) into the thread of the screw tube.

### Please note:

No glue on the outer side of the screw tube is allowed.

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Step 2:

- Move temple in closed position into the slot of the lug.

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Please note:

The of groove the temple hinge is positioned in front direction, that the temple can slide over the threaded pin.

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Step 3:

- Mount screw tube with hinge washer and turn it clockwise with screwdriver until screw head touches the temple.
- Wait until glue is hand tight hardened (hardening time depending on used glue) and then adjust hinge friction.

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Please note:

The screw must not be tightened too strong. Threaded pin can be damaged.

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## 5.7 PORSCHE DESIGN ABSORBER

### Shockabsorber Change

The technical notes refer to the following spectacle frames:

P8973, P8974, P8765, P8766, P8767

For shockabsorber change a thin pointed object is required.

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### Shockabsorber dismounting

Step 1:  
Move temple to closed position

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Step 2:

Use a sharp object to remove the shock absorber from the temple. Remove the pin.

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Step 3:

In case a piece breaks off at the end of the shock absorber and gets stuck in the temple, you can loosen the rest with the help of a needle at the opening provided for this purpose.

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## Shockabsorber mounting

### Step 1:

Put the new shock absorber back into the temple with the peak forward.

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### Step 2:

Close the temple. When you hear a small click, the shock absorber is snapped into its groove.

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## 6.1 FUNCTION

### **UV Radiation**

At high intensity, UV radiation can lead to conjunctivitis and other damage to the eye. Therefore, all Rodenstock sun protection lenses provide a broad protection against UV radiation.

### **Glare Protection**

Glare is reduced by filtering certain colour ranges of the incident light – to varying degrees depending on the tint of the lenses. The percentage value (e.g. 85%) behind the tint indicates the proportion of incident light which is absorbed by the lenses and therefore does not reach the eye.

### **Contrast Enhancing**

In bright sunlight, a high proportion of blue scattered light overlays the other colours and reduces their contrasting effect. Contrast enhancing sun protection lenses filter or minimize the high energy, short-wave blue light from the incoming light and therefore ensure a high-contrast visual experience.

## 6.2 LENSES

### **LAMBDA LENS TECHNOLOGY**

Rodenstock sunglasses with Lambda Lens Technology protect the eyes and optimize colour perception for optimal vision. By filtering disturbing blue stray light, the high-precision sunglass lenses increase contrast and the wearer benefits from greater visual acuity and detail sharpness.

### **Polarized lenses**

With polarized lenses glare is eliminated which leads to a visual clarity and a broad UV protection is guaranteed.

Since reflections are reduced for example on a wet street, a higher safety is given when driving. Improved visual comfort is guaranteed and a better contrast is given which allows a true perception of colours and details. Polarized lenses are ideal for all kinds of sports especially for water sports for example.

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## **VISION DRIVE LENS TECHNOLOGY™**

Thanks to their innovative VISION DRIVE™ lens technology, Porsche Design sunglasses are the perfect companions for everyday life. Especially behind the wheel, the lens technology gives you decisive optical advantages as well as the ultimate driving experience. Getting the best possible view from the cockpit is an important safety factor. Thanks to their excellent quality and image properties, you maintain vision on the road. Even in the most difficult weather and road conditions.

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## **6.3 MATERIALS**

### **Polycarbonate**

The use of polycarbonate for lenses makes them ideal for all kinds of activities due to the extreme resistance to breakage. Polycarbonate offers ideal protection for sports such as squash, tennis or other outdoor activities. Polycarbonate is also highly resistant to blows and impacts.

### **CR-39**

CR-39 is identical with the index 1.5 material and has a small index of refraction. CR-39 lenses are immune to heat and can easily be coloured in all possible shades, which is why they are very often used for sunglasses. Edging CR-39 lenses is very easy to carry out.

### **Polyamide**

Lenses made from polyamide offer real comfort due to their light weight. They are also very robust and highly resistant to breakage, which ensures that they have an excellent service life. The flexibility and versatility is multi-faceted, as the material can be used with any frame material (acetate, metals or injection moulding).

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
## 7.1 TEMPLE TIPS (HEAT-SHRINK TUBING)

### Step 1:

- Remove the temple tip by cutting the heat-shrink tubing lengthwise with a sharp cutter or scalpel
- Pull off heat-shrink tubing

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### Step 2:

- Centering of the new heat-shrink tubing

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
### Step 3:

Shrink the heat-shrink tubing using a frame heater:

- Hold the heat-shrink tubing so that it cannot slip out of place
- Heat the entire heat-shrink tubing even until the tube diameter is halved
- Then shrink the two ends completely and finally the middle part

**Attention:** During the whole process, the temple tip must always be moved and turned over the frame heater to avoid too high temperatures in one place! Touch the heated parts only with suitable protective gloves!

- Finally, allow the temple tip to cool down well

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
## 7.2 LENS FIXATION WITH RIVETS


Required equipment:

- Screw tweezers
- Ball cutter (Ø 0.8-1.0 mm)
- Soft base (e.g. soft tissue)
- Clavulus (Eyeglass rivet press) with rivet head maker
- Spare lens and rivet

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### Step 1:

- Place ball cutter centered on the rivet
- Mill into the rivet and remove it
- Remove the old lens

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
### Step 2:

- Inserting the new lens
- Insert the new rivet into the hole with tweezers

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### Step 3:

- Place the soft base on the anvil of the Clavulus
- Position the frame on it
- Riveting with the Clavulus

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Rodenstock: SERRE001

Porsche Design: SERPD001

The product range is rounded off by additional equipment to make working with Rodenstock and Porsche Design frames considerably easier.

### **Service-Box**

We offer boxes with spare parts. These repair and replacement sets contain a selection of screws, temple tips and nose bridges for the current spectacle frames of Rodenstock and Porsche Design. We offer spare parts for our spectacles for at least 2 years after the date of purchase and frequently even longer – significantly more than what the competition offers.

Rodenstock: SERRE001

Porsche Design: SERPD001

DE RIGO VISION S.p.A.  
Zona Industriale Villanova 12  
32013 Longarone (BL)  
Italy  
[www.derigo.com](http://www.derigo.com)